- de** Gebrauchsanweisung
- en** User's Manual
- fr** Mode d'emploi
- it** Istruzioni per l'uso
- es** Instrucciones de empleo
- pt** Manual de operação
- ru** Руководство по применению


SwiSto3

de	SwiSto3 Iontophoresegerät	4
en	SwiSto3 Iontophoresis device	16
fr	SwiSto3 Appareil d'iontophorèse	28
it	SwiSto3 Apparecchio per iontoforesi	40
es	SwiSto3 Aparato de iontoforesis	52
pt	SwiSto3 Dispositivo de iontoforese	64
ru	SwiSto3 Прибор для проведения ионтофореза	76

de

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


User's Manual

SwiSto3 iontophoresis device for home and clinical use

Dear Customer, thank you for choosing a KaWe product. Our products are known for their high quality and longevity. This KaWe product fulfils the requirements of EG Directive 93/42/EWG (standard for medical products). We assume no liability for improper or inadequate use of this device, or if the device is not checked before use. Please read these instructions carefully and completely before use and follow the maintenance instructions.

1. Application

The SwiSto3 is used in medical facilities and in the home by healthcare professionals as well as by non-professionals.

2. Purpose

The SwiSto3 iontophoresis set is designed to treat intense foot, hand and underarm perspiration. This is a proven and tested anti-perspiration system that uses tap water iontophoresis. A reduction of perspiration can be observed after as few as approx. 15 – 20 treatments, lasting 20 to 30 minutes each.

Table of Contents

	Page
1. Application	17
2. Purpose	17
3. Scope of delivery: SwiSto3 – iontophoresis set	17
4. Special accessories	17
5. Manufacturer responsibility	17
6. Graphical symbols and warnings	17
7. Contraindications / side effects and preventative measures	18
8. Start-up / preparation	19
9. Treatment of foot or hand perspiration	19
10. Treatment of underarm perspiration	20
11. Important information	21
12. Troubleshooting	21
12.1 Treatment does not start – possible causes	21
12.2 Treatment is interrupted – possible causes	21
12.3 Treatment current feels too weak – possible causes	21
12.4 Display does not illuminate – possible causes	21
13. Recurring safety inspections pursuant to § 6 MP regulations for SwiSto3 devices	21
14. Cleaning and disinfection	22
15. Storage while not in use	22
16. Transport and storage	22
17. Disposal	22
18. Technical information	23
19. Guarantee	23
20. Guidance and manufacturer's declaration – electromagnetic immunity	23
21. Recommended separation distances between portable and mobile RF communications equipment and the SwiSto3	26

3. Scope of delivery: SwiSto3 iontophoresis set

- 1 x SwiSto3 iontophoresis device
- 2 x water trays
- 2 x plate electrodes (150 x 220 mm) incl. cable
- 2 x electrode sponge pads
- 1 x charging plug

4. Special accessories

- 2 x plate electrodes with sponge pockets for underarm perspiration treatment, 2 x replacement sponge pockets (ca. 90 x 110 mm) alone

5. Manufacturer responsibility


The manufacturer only assumes responsibility for the safety, reliability and efficient operation of the SwiSto3 iontophoresis device if:

1. repairs, modifications, readjustments or enhancements are performed by the manufacturer or by a person authorised by the manufacturer and
2. the electrical system installed in the treatment room meets the requirements of ICE-stipulations and
3. the SwiSto3 iontophoresis device is used in accordance with the instructions given in the User's Manual.

No liability will be assumed for damages resulting from improper use.

Only use the plate electrodes and connection cables supplied by KaWe and/or the original KaWe-brand spare parts and accessories.

6. Symbol key:

	Power button
	Anode
	Cathode
	Increase current
	Decrease current
	Menu button for selecting settings
	Continuous current
	Pulsed current
	Treatment time
	Battery symbol
	Electromagnetic field warning
	Charge symbol
	BF-type device according to DIN IEC 601 Part 1/VDE 0750 Part 1
	Heed the User's Manual
	Caution!

	Separate disposal of electric and electronic devices
	Manufacturer
	Date of manufacture
	Lot code
	Complies with relevant EU guidelines
iP41	Protected against water droplets and solid matter with a diameter of greater than or equal to 1.0mm
	Product reference number
	Series number
	Temperature limit
	GOST-R certification of exports to Russia

7. Contraindications/side effects and preventative measures

IMPORTANT! Vital information:

- Consult your physician before using this device to receive an individualised treatment plan. Tap water iontophoresis may not be performed under the following circumstances or if the patient has one of the following conditions (contraindications):

- implanted electronic devices (i.e. pacemaker)
- pregnancy
- metal implants in the areas of the body through which the current is conducted (arms or legs)
- metal-containing intrauterine devices (contraceptive coils) when treating the feet
- large-sized skin defects that cannot be covered with Vaseline or isolating bandages/wraps
- heart murmurs or insensitivity to pain.

The device may only be installed and operated in accordance with the instructions given in the user's manual. Radios, mobile telephones or other similar devices that may affect the device must be placed at least 2 metres away from the device. The device may not be used by unsupervised children and must be stored out of reach of children.

The cables pose a risk of strangulation. Small parts may become loose and could be swallowed or inhaled.

- Only use standard accessories provided by KaWe. The use of other cables, charging devices, electrodes, water trays etc. is not allowed. The use of non-standard parts can influence the EMC interference immunity and electrical safety of the unit.
- Simultaneous treatment with another high-frequency surgical instrument may cause skin burns under or near the electrodes.
- The operation of the iontophoresis device within close range (approx. 1 m) of a short-wave or a microwave therapy unit can cause fluctuations in the output readings of the iontophoresis device and therefore, should be avoided.

- Never perform treatments during thunderstorms. If necessary, immediately discontinue the treatment and turn off the unit.
- The maximum current and the maximum voltage are limited by safety values that are determined by the regulations for electrical medical devices, thereby eliminating the risk of harm to the patient.
- Electrical shocks are primarily caused by contact failures between the connector cables and the electrodes. Over time, transfer resistance occurs between the electrodes and the connector cables due to the direct contact with water. This can be avoided by disconnecting the cable on the electrode just before beginning treatment and then reinserting it with a twisting motion. The friction thereby produced creates a clean metal-to-metal contact which is essential for a constant flow of current. This practically eliminates uncomfortable current peaks that could otherwise be felt by the patient.
- Persons with prosthetics must remove these before treatment. They must also consult a doctor before using this device.
- It is essential that the patient ensure that the success of the treatment is not affected by any type of disturbance. Due to this reason as well for electrical safety reasons, it must be ensured that no children or house pets are allowed in the room at the time of treatment.
- The patient should be aware of the fact that – especially during underarm perspiration treatment – any rash movement during the treatment can alter the surface on the body being used for the current transfer. As a result, the current load on the affected body part may exceed the allowed value, which although it cannot cause actual harm, could potentially cause discomfort, skin irritation and/or slight burns. Such complications can be best avoided by:
 - Avoiding direct metal contact between the electrodes and the body by using the foam pads or sponge pockets.
 - Removing the hands and/or feet slowly from the water trays while the treatment current is flowing.
 - Ensuring that the electrodes are always fully inserted into the sponge pockets during underarm perspiration treatments.
 - Ensuring that during underarm perspiration treatment the contact pressure between the sponge pocket and the underarm remains as constant as possible for the entire treatment time.
- Direct current tap water iontophoresis is not suitable for sensitive users as it may cause skin irritation or paresthesia. In order to minimize these side effects, this device is to be used with pulsed current and the current must not exceed given amperage. Nevertheless, while using pulsed current, slight stinging, burning or tingling sensations or a slight reddening of the treated skin areas may still occur. In order to prevent these side effects, the sponge pockets should only be lightly and evenly moistened and not soaked before they are placed under the armpits.
- Cover damaged/rough areas of the skin on the palms of the hands, the soles of the feet and the grooves of the nails with Vaseline or another greasy ointment, since the electrical current transmissibility through these areas is increased.

- Using this device too often may lead to minor skin irritation in some cases.
- If the treatment is intentionally performed with a very high electric current, an uncomfortable tingling can occur in the extremities through which the current is flowing.
- The complete surface of the body parts being treated must rest with uniformly-distributed pressure on the foam pads and/or sponge pockets. A non-uniform distribution of the surface and of the pressure may cause a current density ($>0.2\text{mA}/\text{cm}^2$) that is too high in some areas, which can cause skin irritation or minor burns.
- Jewellery and wrist watches should be removed before treatment.
- The optimum therapy results are achieved by selecting the highest possible (however, different for each individual) current intensity in the water bath during the treatment.
- The influence on sweat inhibition is much more pronounced on the positive (+) pole than on the negative (-) pole. Therefore, if applicable, it is recommended that the patient change the polarity of the electrodes after each treatment.
- The hands and the feet cannot be treated simultaneously. In accordance with the user's manual, hand, feet and underarm treatments are to be carried out as separate procedures.
- In order to avoid skin irritation or burns, only electrodes recommended by KaWe are to be used.
- Please ensure that the included foam pads are always placed on the electrode plates!
- In order to reduce the skin's natural resistance and if necessary, to make the treatment current feel more comfortable to the user, the hands or the feet should be massaged in a warm water bath prior to being treated.
- Renew the water in the trays before each use.
- Ensure that the display can be read during treatment and is not illegible due to lighting conditions.
- The device and its parts (such as the membrane keypad etc.) are to be checked before each use to ensure that they are in sound condition. Should there be any apparent problems such as wear etc., please send the device in to the manufacturer.
- The expected operational life time of the device is, depending on use and maintenance, between 4 and 8 years. The built-in rechargeable battery is designed for 800 complete charging cycles. With a fully-charged battery, approx. 10, 20 minute-long treatments can be performed at full power.
- If necessary, non-expert users are to contact the manufacturer or their representative in order to receive help with the device or to report unexpected functional operation or occurrences.

8. Start-up/preparation:

- Caution, condensation can cause a safety hazard! If the device becomes very cold during storage or transport, ensure that it reaches room temperature before using it.
- When operating the device, heed the instructions given in "Guide-

- lines and manufacturer declaration – Electromagnetic immunity".
- To charge the battery, first plug the low-voltage connector into the socket on the back side of the unplugged SwiSto3 SwiSto3 iontophoresis device. The plug-in charger is to be stored in a location that is easily accessible and easy to reach.
- Ensure that the connector cable between the adapter and the iontophoresis device is not a tripping hazard!
- Charge only with the KaWe REF 05.19170.002 charging plug!
- Charge the battery completely before using for the first time. The charging time for a completely discharged battery is approx. 10 hours. When charging, the battery symbol blinks and the charging symbol on the display appears in the middle of the screen. The device cannot be used while it is being charged. The charging progress can be seen on the battery symbol.
- To perform a treatment, place the SwiSto3 iontophoresis device on a flat and dry surface.
- Place the water trays in an appropriate location.
- Place the electrodes (REF 05.19040.021) in the water trays.
- Place the foam pads on the electrodes such that the electrodes are completely covered.
- Use the connector cables to connect the electrodes with the output jacks (+ and -) on the side of the device. When connecting the cables to the electrodes, they should be plugged in with a simultaneous twisting motion.
- Fill both of the tubs with a minimum of 3-4 cm of lukewarm tap water. Take special care that the trays do not overflow or have the potential to overflow when the hands or feet are put into the water.
- For the underarm perspiration treatment, insert the plate electrodes into the sponge pockets and soak them well in water.
- Power-on the SwiSto3 iontophoresis device by pressing the push-button switch. Hold the button down for approx. 1 second until the display illuminates. Current treatment values will be shown on the display.
- The SET button is used to select each of the individual settings. Each of the blinking values can be changed with the up and down buttons. The following can be changed: the current intensity in steps of 1 mA, continuous current / pulsing current and the treatment time in min:sec.
- If the display is not readable or the unit appears to be damaged in any way, send the device in to the manufacturer immediately.

9. Treatment of foot or hand perspiration

- Prepare device for use as described in Section 8 "Start-up/preparation".
- Power-on the device by pressing the power button  for one second.
- In each of the tap water-filled trays, place the palm of a hand or the sole of a foot onto the foam pad-covered plate electrode. Directly touching the electrodes is not dangerous; in order to avoid spots of excessive current density, it is recommended however, that as

much of the surface of the body part to be treated as possible rests on the electrodes and the pressure is uniformly distributed. The parameters that can be selected blink on the display. The current intensity for the treatment can be selected with the UP and DOWN buttons. --> Confirm selection by pressing the SET button. --> Select pulsing current with the UP button and continuous current with the DOWN button. --> Confirm selection by pressing the SET button. --> Set the treatment time in minutes (MIN) using the UP and DOWN buttons. --> Confirm selection by pressing the SET button. --> Set the treatment time in seconds (SEC) using the UP and DOWN buttons. --> Confirm selection by pressing the SET button. --> Setting selection is completed when the treatment starts, the time countdown begins and the clock symbol begins to blink. To begin the treatment, complete the electric circuit by putting your hands or feet into the water trays.

- By pressing the  button on the membrane keypad on the upper side of the device the treatment current can be increased until a mild to strong tingling can be felt. The intensity of the current can be seen on the display. Similarly, the current can also be reduced during the treatment by pressing the  button. When the current is being changed during the treatment, the current reading on the display does not blink.
- The treatment stops automatically after the set treatment time has been reached or when the current circuit is broken before the set time has elapsed. The treatment will resume for the remainder of the set treatment time when the current circuit has once again been closed. The clock symbol  blinks and the time countdown resumes.
- If the clock showing the treatment time  has run down to 00:00 and the current circuit is still complete, the + / - symbols on the display both begin to blink. Once the set treatment time has expired, none of the settings can be changed until the circuit has been interrupted by removing the hands or feet from the water trays. The + / - symbols then stop blinking and the MIN symbol begins to blink, signaling that new settings can now be selected. The previously set treatment time will be saved.
- Turn the SwiSto3 iontophoresis device off after completing treatment by pressing the power button . The device will also turn off automatically after two minutes.
- The last-used treatment settings are automatically saved and used when the device is restarted for the next treatment. These settings can be changed with an open current circuit as described in Section 9.
- The SwiSto3 iontophoresis device can also be plugged into the local power supply network to be charged using the charging plug. The battery symbol  shows the battery charge up in the top right-hand corner of the display.
- The charging symbol  appears when the charging plug is plugged in. After testing the rechargeable battery for at most one minute, the charging function begins and the battery symbol begins to blink.
- When the charging plug is plugged in, the treatment is interrupted.

10. Treatment of underarm perspiration

- Prepare device for use as described in Section 8 "Start-up/preparation".
- Power-on the device by pressing the power button  for one second.
- For the underarm perspiration treatment, **wet the sponge pockets evenly with tap water.** In order to avoid the occurrence of current densities that are too high on isolated spots on the body however, ensure that as large a contact surface as possible is touching the electrodes with a pressure distribution that is as uniform as possible.
- Insert the plate electrodes (REF 05.19080.021) into the sponge pockets.
- The parameters that can be selected blink on the display. The current intensity for the treatment can be selected with the UP and DOWN buttons. --> Confirm selection by pressing the SET button. --> Select pulsing current with the UP button and continuous current with the DOWN button. --> Confirm selection by pressing the SET button. --> Set the treatment time in minutes (MIN) using the UP and DOWN buttons. --> Confirm selection by pressing the SET button. --> Set the treatment time in seconds (SEC) using the UP and DOWN buttons. --> Confirm selection by pressing the SET button. --> Setting selection is completed when the treatment starts, the time countdown begins and the clock symbol begins to blink. To begin the treatment, complete the electric circuit by placing the sponge pockets with the plate electrodes in them under your armpits.
- Raise the upper arm and place the sponge pocket on the armpit.
- The sponge pocket is then held in place by pressing the upper arm down against the body. Before lowering your upper arm ensure that the electrodes and connector cables are inserted completely into the sponge pockets, as direct contact between the metal and the skin could cause localised burns. Especially when performing an underarm perspiration treatment, ensure that the pressure exerted by the upper arm remains as constant as possible for the entire duration of the treatment – which can be up to 30 minutes.
- By pressing the  button on the membrane keypad on the upper side of the device the treatment current can be increased until a mild to strong tingling can be felt. The intensity of the current can be seen on the display. Similarly, the current can also be reduced during the treatment by pressing the  button. When the current is being changed during the treatment, the current reading on the display does not blink. **WARNING:** Do not exceed the recommended maximum current intensity of 5 mA.
- The treatment stops automatically after the set treatment time  has been reached 00:00 or when the circuit is broken before the set time has elapsed. The treatment will resume for the remainder of the set treatment time when the circuit has once again been closed.
- If the clock showing the treatment time  has run down to 00:00

and the current circuit is still complete, the + / - symbols on the display both begin to blink. Once the set treatment time has expired, none of the settings can be changed until the circuit has been interrupted by removing the electrodes. The + / - symbols then stop blinking and the MIN symbol begins to blink, signaling that new settings can now be selected. The previously set treatment time will be saved.

- Turn the SwiSto3 iontophoresis device off after completing treatment by pressing the power button . The device will also turn off automatically after two minutes.
- The last-used treatment settings are automatically saved and used when the device is restarted for the next treatment. These settings can be changed with an open current circuit as described in Section 9.
- The SwiSto3 iontophoresis device can also be plugged into the local power supply network to be charged using the charging plug. The battery symbol shows the battery charge up in the top right-hand corner of the display.
- The charging symbol appears when the charging plug is plugged in. After testing the rechargeable battery for at most one minute, the charging function begins and the battery symbol begins to blink.
- When the charging plug is plugged in, the treatment is interrupted.

11. Important information

If the treatment circuit is broken, the current is automatically shut off so that electric shock does not occur. When the battery charge is low, the battery symbol begins to blink. Treatment is only to be carried out after first charging the battery. When the battery is empty, the device turns off automatically.

If the rechargeable battery has an extremely low charge, the display blinks while charging. After being charged for max. 30 minutes, the display will illuminate normally without blinking.

Recommended maximum values	mA
Underarm perspiration treatment	5
Hand perspiration treatment	15
Foot perspiration treatment	25

12. Troubleshooting

12.1 Treatment does not start – possible causes

- Treatment circuit is not closed. Eliminate possible contact failure between the connector cables and the electrodes by unplugging the electrode connection pin and reinserting it into its socket with a twisting motion. Then close the treatment circuit through your body by placing your hands or feet into the water onto the foam pad-covered electrodes, or when treating underarm perspiration, by placing

the sponge pockets under your arms.

- Sponge pockets (for underarm perspiration treatment) are not wet enough. The entire surface of the sponge pockets must be uniformly moistened.
- Too little tap water in the trays. The level of the water in the trays must be approx. 3-4 cm.
- The electric conductivity of the tap water is not sufficient. In this case, the conductivity of the tap water must be increased by adding minerals to the water.
- The battery charge is low. The battery symbol blinks, + and - symbols blink when the treatment circuit is complete.

12.2 Treatment is interrupted – possible causes

- Disruption of the treatment circuit – such as by suddenly removing hands/feet from the trays or raising your arm during an underarm perspiration treatment.
- Electrode contact failure. Eliminate possible contact failure between the connector cables and the electrodes by unplugging the electrode connection pin and reinserting it into its socket with a twisting motion.
- The battery charge is low. The treatment current is turned off. The battery symbol blinks, + and - symbols blink when the treatment circuit is complete. Charge the battery as described in section 8.

12.3 Treatment current feels too weak – possible causes

- Electrodes are worn. After prolonged use, the condition of the electrodes can deteriorate (due to calcium deposits on their surface, for example) to such an extent as to inhibit the power output of the device. The electrodes are to be replaced.

12.4 Display does not illuminate – possible causes

- Battery is completely discharged. First charge the battery with the supplied charging plug. The charging symbol will appear on the display. Soon after the battery has begun charging, the device switches automatically to charging mode.
- If the charging plug is plugged in and still no charging symbol appears on the display, check the voltage of the electrical outlet that is being used with appropriate equipment. If the electrical outlet is functioning properly, the SwiSto3 and the charging plug must be inspected at a location authorised by KaWe – see as well the Section “Safety Inspections”.

13. Recurring safety inspections pursuant to § 6 MPG regulations for SwiSto3 devices

For the SwiSto3, the following inspections must be performed at least every 12 months by individuals able to perform such safety inspections properly due to their training, knowledge and practical experience, and do not require any special instructions with regard to these

inspection activities.

- Visually inspect the device and accessories for mechanical damage that would adversely affect their functionality.
- Check safety-related labels and the display to insure that they are legible.
- Perform function check as instructed in the user's manual.
- Test the iontophoresis current with a load resistance of 1.5 K Ω .
- The SwiSto3 device contains no parts that require maintenance by the user.
- The housing may not be opened.

When the UP button  is pressed, the current intensity is set to 25 mA. For this, a current of 25 mA \pm 2 mA is necessary. Then a short-circuit test is performed on the output, whereby the values must not change.

When the DOWN button  is pressed, the current intensity output of the device gradually diminishes. Here, the current value must decrease with each decreasing display increment. The treatment current is tested with a reference instrument.

During technical safety inspections, an electrical safety inspection pursuant to the German standard DIN VDE 0751-1; VDE 0751 Part 1 (Repeat Examinations and Tests Prior to Operation of Medical Electrical Equipment or Systems - Part 1: General Requirements) or DIN EN IEC 60601-1 is to be performed.

The equivalent leakage currents may be up to 1.5 times that of the first measured value and at the same time may not exceed the limiting value. The first measured values are given in the enclosed test reports. The safety inspection must be entered into the technical manual in accordance with § 6 MPG of the regulations pertaining to the sale of iontophoresis devices and the inspection results must be documented. If the device is found to not be functional and/or safe to operate, it is to be repaired or the user of the device is to be notified of the risks of using the device.

14. Cleaning and disinfection

- Empty the water trays.
- Wipe the water trays, electrodes and connector cables with a dry cloth and then disinfect them as follows: use any commonly-used disinfecting agent such as Bode Bacillol[®], Braun Meliseptol[®] or Orochemie B30. The instructions given by the manufacturer of the disinfecting agent are to be followed. The surface that is to be disinfected must be wiped with the disinfectant using light pressure, such that the surface is wet enough and enough of the active ingredient of the disinfecting agent is deposited. The applied disinfectant solution is to be allowed to dry on the surfaces and they may not be

wiped dry. The reaction/wait time for the above-named disinfecting agents is 5 minutes and must be observed

- The SwiSto3 Iontophoresis Device is to be wiped down weekly with a disinfecting agent.
- Water trays, electrodes and connector cables are to be disinfected before each use! (Wipe or spray disinfectant)
- Sponge pockets can be rinsed with warm water, but not with any kind of cleaning agent and for hygienic reasons, are each to be used only by one and the same patient.

Important!

Parts such as electrodes, foam pads, sponge pockets and connector cables are subject to a certain amount of wear and may require replacement after prolonged use. Check all parts for damage before each application, especially the connecting lines. If corroded, the cables are to be replaced. If the contact plugs on the SwiSto3 become corroded, send the device in for inspection.

The integrated rechargeable battery may only be replaced by the manufacturer. Should the battery be defective (will not charge, discharges too quickly), send the device in for inspection/repair.

Warning: Modifying the SwiSto3 ME device is not allowed! In order avoid the risk of electric shock, only the included charging plug (REF 05.19170.002) may be used to charge the unit.

15. Storage while not in use

In order to prevent the SwiSto3 iontophoresis device from becoming damaged, it should be stored along with its accessories in its original packaging if not in use for a prolonged period of time.

16. Transport and storage

When transported, the SwiSto3 iontophoresis device is to be packaged such that it is protected from damage. If the device is not to be used for a prolonged period of time, it should be charged completely before being stored. The device and accessories are to be stored in a dry location between 0 °C and +40 °C.

17. Disposal

The non-accessible part of the SwiSto3 iontophoresis device contains electrical and electronic components. For this reason, at the end of its service life, the device must be brought to an appropriate waste disposal site in accordance with local regulations or returned to the manufacturer.

18. Technical Information

Classification	Applied part: type BF
Protection classification	IP 41
Dimensions	W 110 x H 84 x L 120 mm
Nominal output current	Adjustable 1 25 mA
Load resistance	max. 1.5 kΩ
Power input	1.5 W
Power supply	Li-ion rechargeable battery
Operation mode	Continuous operation
Weight	Charging plug: approx. 0.10 kg iontophoresis device: approx. 0.36 kg
Safety features	Electronic current limiting device for the treatment current, Current is cut off when the treatment circuit is broken, System interrupt for when the treatment circuit is not closed
Current type	Continual current or pulsed current up to a max. of 25 mA

Ambient conditions

Operating temperature	+ 10 °C to + 40 °C
Storage temperature	0 °C to + 40 °C
Relative humidity	30 % to 75 %
Air pressure	700 hPa to 1060 hPa

20. Guidance and manufacturer's declaration – electromagnetic emissions

The SwiSto3 is designed to be used in the types of electromagnetic environments listed below. The customer or user is responsible for ensuring that this device is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – Guidance
RF emissions CISPR 11	Group 1	The SwiSto3 device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The SwiSto3 device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

19. Guarantee

We guarantee this product (except for the Li-ion rechargeable battery) for 48 months. The prerequisite for this guarantee is a proper use of the product in accordance with this instruction manual. The guarantee is valid from date of purchase. The bill of sale serves as proof in this case. Wear parts – such as electrodes, foam pads and sponge pockets – are not covered by the warranty!

Before returning the instrument in the event of a malfunction, the possible causes described in the section “Troubleshooting” should all be examined thoroughly. Should we discover upon inspecting the returned device, a reason for the malfunction that is described in the user's manual that could have been easily solved by the customer themselves, we will charge an inspection fee. (This fee includes the cost of labour, postage and packing, as well as VAT.)

EMC instructions: Medical electrical equipment is subject to special safety measures so as to ensure safe operation with respect to electromagnetic compatibility. The following information helps to ensure safe operation.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The SwiSto3 is intended for use in the electromagnetic environment specified below.
The customer or user is to assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1kV voltage line to line ± 2kV voltage line to ground	± 1kV voltage line to line ± 2kV voltage line to ground	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power-supply input lines IEC 61000-4-11	<5 % U_T for 1/2 cycle (>95 % dip) 40 % U_T for 5 cycles (60 % dip) 70 % U_T for 25 cycles (30 % dip) <5 % U_T for 5 seconds (>95 % dip)	<5 % U_T for 1/2 cycle (>95 % dip) 40 % U_T for 5 cycles (60 % dip) 70 % U_T for 25 cycles (30 % dip) <5 % U_T for 5 seconds (>95 % dip)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued OPERATION during a power mains interruption, it is recommended that the device be powered by an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	If interference occurs, it may be necessary to place the SwiSto3 further away from the sources of power frequency magnetic fields or to install a magnetic shield: The power frequency magnetic field should be measured at the designated site to ensure that it is not large enough to cause interference.

Note: U_T is the ac mains voltage prior to application of the test level.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The SwiSto3 is intended for use in the electromagnetic environment specified below. The customer or user of the SwiSto3 is to assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the SwiSto3 device, including cables, than the recommended separation distance, which is calculated using the equation for the frequency transmission.</p> <p>Recommended separation distance $d = 1.2\sqrt{P}$</p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	<p>$d = 1.2\sqrt{P}$ for 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ for 800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (M).</p> <p>Field strengths^a from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level^b in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a: Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the SwiSto3 device is used exceeds the applicable RF compliance level above, the SwiSto3 device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

21. Recommended separation distances between portable and mobile RF communications equipment and the SwiSto3

The SwiSto3 is intended for use in the electromagnetic environment described below. The customer or the user of the SwiSto3 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the SwiSto3 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance (meters) according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.10	0.38	0.38	0.73
1	1.20	1.20	2.30
10	3.80	3.80	7.30
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (w) according to the transmitter manufacturer.

Note 1: For the calculation of the separation distance for transmitters with a frequency range of 80 MHz to 2.5 GHz, an additional factor of 10/3 is used in order to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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KIRCHNER & WILHELM GmbH + Co. KG